







	JOBST FarrowWrap 4000	JOBST FarrowWrap Lite	JOBST FarrowWrap Strong	JOBST FarrowWrap Classic
Style				
Compression Level (mmHg)	30-40	20-30	30-40	30-40
Indications	For venous disease with or without leg ulcers	Soft fabric with reduced compression , suitable for palliative patients or those with lower ABPI's	Soft fabric with strong compression , suitable for patients with some shape distortion	Rigid fabric for patients with more severe oedema and irregular limb shape, fabric lays flat against folds
Level of Oedema	Mild	Mild to Moderate	Mild to Severe	Mild to Severe
Sizes	Extra Small, Small, Medium, Large	Extra Small, Small, Medium, Large, Extra Large		
Length	Regular (30-33cm) Tall (34-37cm)	Foot piece = Regular & Tall Leg piece = Regular & Tall Thigh Piece = Short, Regular, Tall		
Bands	4 x interlocking bands with integral sleeve	50% overlapping bands to imitate short stretch bandaging		
Liners	1 x JOBST FarrowHybrid ADI. (Available to prescribe individually in 20-30 (ADI) or 30-40 (ADII) mm/Hg in the foot only)	Beige Liner (Small or Large) provided with foot and leg pieces TG Soft Liner (Small or Large) provided with thigh piece Provides a smooth base/skin protection layer when applying the outer wrap pieces		