

Tissue Viability Clinical Decision-Making Tool for Compression Hosiery for **Normal Limb Shape** (Ready-To-Wear / Off-The-Shelf)




Level of oedema	None		Mild		Mild to Moderate			Moderate
Compression Level mmHg	10mmHg	British Standard CCL1 - 14-17mmHg CCL2 - 18-24mmHg CCL3 - 25-35mmHg	40mmHg (Combined) 1 x liner 1 x stocking	RAL CCL 1 - 18-21mmHg CCL 2 - 23-32mmHg	RAL CCL 1 - 18-21mmHg CCL 2 - 23-32mmHg			RAL CCL 2 - 23-32mmHg CCL 3 - 34-46mmHg
Ready to Wear (Circular knit)	Altipress Liners	DuoMed Soft	Altipress 40	JOBST UltraSheer	JOBST Opaque	JOBST forMen	Mediven for Men	JOBST Classic
Garment information	Palliative care, maintenance or those intolerant to higher levels of compression	Daily wear for gentle care	For the treatment of venous leg ulcers	Silky look with sheerness to encourage patient compliance	Soft fitting and high comfort with superior breathability & moisture management	High-quality cotton content for ease of donning and increased comfort	Fibres transports moisture to the outer surface of the sock to maintain comfort	Firmer fabric & stiff resistance with reinforced toe & robust heel providing high durability with a comfortable fit
Colours	White & Beige	Black & Sand	10mmHg White Liner 25-35mmHg Beige Stocking	Natural, Black, Caramel, Bronze, Espresso, Anthracite	Natural, Black, Caramel, Bronze, Navy	Black, Navy, Dark Grey, Brown, Khaki	Black, Navy, Brown, Anthracite	Sahara & Pearl
Styles	Knee High	Knee high & Thigh High	Knee High	Knee high, Thigh High, Tights		Knee High		Knee high, Thigh High, Tights
Lengths	Short, Regular & Long	Petite & Long	Short, Regular & Long	Petite & Regular		Regular & Long		Petite & Regular
Toe options	Closed	Open & Closed	Closed Toe Liner Open Toe Stocking	Open & Closed		Closed		Open & Closed

For advice on hosiery please complete a lower limb referral and send to gwh.schs.tvteam@nhs.net

Alternatively, contact Dan Young at Essity on Daniel.young2@essity.com / 07711466482 or Jo Driscoll at Medi on jo.driscoll@mediuk.co.uk / 07795668403

or Katie Rose-Burrows at Urgo on k.rose-burrows@uk.urgo.com / 07966229427

Tissue Viability Clinical Decision-Making Tool for Compression Hosiery
for **Abnormal Limb Shape** (Made to Measure / Custom Fit)

Level of Oedema	No Oedema	Mild to Moderate	Any
			
Compression Level	10mmHg	RAL CCL1 - 18-21mmHg CCL2 - 23-32mmHg CCL3 - 34-46mmHg	RAL CCL1 - 18-21mmHg CCL2 - 23-32mmHg CCL3 - 34-46mmHg
Made to Measure (Flat Knit)	Altipress Liners (Circular knit)	JOBST Elvarex Soft	JOBST Elvarex
Indications	Palliative Care	Ideal for patients with dry, sensitive or fragile skin Latex free	Robust material providing high resistance for stubborn & rebound oedema Contains Latex
Colours	White	Beige, Black, Cranberry, Dark Brown, Grey, Dark Blue, Ruby Red	Beige, Black, Cranberry, Dark Brown, Grey, Dark Blue, Henna, Denim, Graphite
Styles	Knee High	Knee High, Thigh High, Tights, Bermuda, Capri Tights	Knee High, Thigh High, Tights, Bermuda, Capri Tights, Chaps
Lengths	Made to Measure	Made to Measure	Made to Measure
Toe Options	Closed	Open & Closed	Open & Closed
Options	-	Silicone Band, SoftFit, T-heel, Silk Pocket	Zipper, Silicone Band, SoftFit, Ankle Pads, T-heel, Silk Pocket

For advice on hosiery please complete a lower limb referral and send to gwh.schs.tvteam@nhs.net

Alternatively, contact Dan Young at Essity on Daniel.young2@essity.com / 07711466482 or Jo Driscoll at Medi on jo.driscoll@mediuk.co.uk / 07795668403

or Katie Rose-Burrows at Urgo on k.rose-burrows@uk.urgo.com / 07966229427